## **Respiratory Virus Panel**

**English**

An indirect immunofluorescence assay for the screening and identification of Adenovirus, Influenza A & B, Parainfluenza type 1, 2 & 3 and RSV.

## Table of Contents

Intended Use	3	Interpretation of Results	22
Introduction	3	Expected Values	24
Assay Principle	7	Limitations of Use	25
Kit Components	8	Performance Characteristics	26
Materials Provided	8	Summary of RVP Procedure	30
Additional Materials Required	10	Interpretation of Symbols	32
Precautions	11	Additional Biotrin Products	33
Safety	11	Bibliography / References	107
Procedural	12		
Storage and Stability	13		
Specimen Collection and Preparation	13		
Assay Procedure	19		

## Intended Use

The Biotrin Respiratory Viral Panel is an in vitro immunofluorescence assay for screening and identification of adenovirus, influenza A, influenza B, parainfluenza 1, parainfluenza 2, parainfluenza 3 and respiratory syncytial virus (RSV) by indirect detection and in cell culture.

## Introduction

Over 200 viruses cause respiratory infection but only seven - adenovirus, influenza A & B, parainfluenza 1, 2 & 3 and respiratory syncytial virus (RSV), are responsible for most of the severe diseases in the very young and immunocompromised host. These viruses are also responsible for significant morbidity in healthy adults during epidemic periods. Nosocomial infections with RSV and influenza, may prove fatal for hospitalised infants and acute care patients.

Early identification of respiratory viruses is essential for effective diagnosis and patient management. The three most common methods of laboratory identification are (a) direct detection, (b) culture confirmation and (c) serology. Of these, direct detection and cell culture confirmation are the standard procedures used by most virology laboratories. The Biotrin Respiratory Virus Panel utilises high quality group-specific and type-specific monoclonal antibodies for direct detection and culture confirmation; providing clear, easy-to-interpret results within 60 minutes.

**Table 1: Seasonality respiratory virus disease**

<b>Virus</b>	<b>Time of Year</b>
<b>Adenovirus</b>	All year
<b>Influenza A</b>	Winter months
<b>Influenza B</b>	Winter months
<b>Parainfluenza 1</b>	All year with peaks in late summer / autumn
<b>Parainfluenza 2</b>	All year with peaks in late summer / autumn
<b>Parainfluenza 3</b>	All year with peaks in late summer / autumn
<b>RSV</b>	Late autumn to early spring

### ***Adenoviruses***

Human Adenoviruses are associated with a wide range of clinical diseases including respiratory, ocular and gastrointestinal tract infections. Infections are particularly common in children and young adults, occurring both sporadically and in outbreaks. In immunocompromised patients severe systemic infections can occur which may be life threatening.

Particular serotypes of adenovirus are associated with the various clinical conditions. Types 15-24 and 37 may lead to ocular disease, ranging from mild conditions such as water borne “swimming pool” conjunctivitis to severe epidemic keratoconjunctivitis. Adenovirus has been associated with 4-15% of all nosocomial cases of viral gastroenteritis in infants; types 40 and 41 are responsible for the majority of these cases.<sup>(1,2)</sup>

Adenovirus is involved in most major types of respiratory disease syndromes. Upper respiratory

diseases caused by adenovirus occur mainly in infants and young children and include colds, pharyngitis and tonsillitis. Adenovirus is also associated with lower respiratory tract infections including bronchitis and bronchiolitis.<sup>(3)</sup> In addition, adenovirus (types 38, 39 and 17) is associated with approximately 5% of acute respiratory disease (ARD) in children and 10% of febrile and childhood pneumonias.<sup>(4)</sup>

Laboratory diagnosis of adenovirus infection plays an important role in effective patient management and the control of outbreaks. Cell culture is the conventional method for the Adenovirus identification and they may be cultured and isolated in a variety of cell lines and identified by a number of methods. Suitable cell lines are Hep-2, HeLa, KB and HEK. The cytopathic effect (CPE) is typically “grape-like” clusters of cells which generally appear within 3 to 10 days. CPE is then confirmed by immunofluorescence, EIA or DNA hybridisation. EM of faeces is also a

standard adenovirus diagnostic method, however this facility is only available in specialised laboratories.<sup>(5,6)</sup>

### ***Influenza A & B***

There are three types of influenza virus; types A, B and C. Influenza A and B are responsible for a highly contagious respiratory disease which occurs during the cold weather months and is often associated with large outbreaks and epidemics. Influenza C rarely results in lower respiratory tract illnesses and is more commonly associated with sporadic upper respiratory tract disease.<sup>(7)</sup>

Influenza virus infection results in abrupt febrile illness and is associated with pharyngitis, laryngitis, croup, bronchitis/ tracheobronchitis as well as bronchitis, influenza and pneumonia. The old and young are partially at risk of severe infection as pulmonary or cardiac complications may result.<sup>(8)</sup>

The classical method for diagnosis of influenza viruses is isolation in PMK, A549 or MDCK cell lines as well as embryonated hen's eggs. CPE appears 3-7 days after inoculation as vacuolation and cell degradation.<sup>(9)</sup>

### ***Parainfluenza 1,2 & 3***

Parainfluenza viruses, combined with RSV are the most important respiratory viral pathogens in infants and children. In older children and adults, illness may be asymptomatic or have symptoms similar to the common cold.

Four types of parainfluenza have been identified. Types 1 and 2 are major causes of croup which is particularly severe in children aged 2-4 years. Parainfluenza 3 may also cause croup but it is more commonly associated with bronchiolitis and pneumonia in infants (most severe infection occurs in children less than 1 year of age). Parainfluenza 4 has been associated with mild upper respiratory tract infections in children and adults.<sup>(7,12)</sup>

Parainfluenza grow well in PMK, LLC-MK2 and HEK cell lines. They are also easily recovered from Vero, A549 and human diploid fibroblasts. CPE appears 4-7 days post inoculation and is observed as an increase in small rounded cells for type 1, syncytial cell formation for type 2 and "stringy edges" for type 3.<sup>(9, 10)</sup>

Confirmation of cell culture has generally been achieved using guinea pig erythrocytes. However, immunofluorescence (IF) with antibodies provide a rapid, inexpensive method capable of type specificity. IF has also been shown to be the best technique for direct detection in respiratory epithelial cells.

### ***Respiratory Syncytial Virus***

RSV is the leading cause of lower respiratory tract infections in infants and younger children. Infection usually results in bronchitis / tracheobronchitis but it may also result in croup, cold and pneumonia. In adults and older children, infection tends to be asymptomatic or expressed

as the "common cold". It is observed in distinct seasonal outbreaks starting in November or December and continuing for approximately 3 months.<sup>(14)</sup>

Direct detection (using IF or EIA) is generally accepted as the method of choice. However, culture confirmation may be used either alone, or as a supplement to direct detection results. For primary virus isolation Hep-2 HeLa are the optimal cell lines, however Vero, LLC-MK2 and CV-1 have also been used. The virus produces a characteristic CPE with syncytium formation and cell destruction. IF is the method of choice as a confirmation immunoreagent.<sup>(15,16)</sup>

### ***Assay Principle***

The **Biotrin Respiratory Virus Panel** utilises an indirect immunofluorescence antibody technique for identifying the 7 major viral respiratory pathogens in infected tissue cultures as well as in direct swab samples. The anti-viral screening reagent is used to confirm the presence of a

respiratory virus. Specific virus identification is then achieved using the anti-adenovirus, anti-influenza A & B, anti-parainfluenza 1, 2 & 3 and anti-RSV monoclonal antibodies.

Fixed patient specimens are incubated with the anti-viral mouse monoclonal antibody on a glass slide. If specific viral antigen, is present a stable complex is formed with the antibody. Following a wash step, FITC labelled goat anti-mouse IgG is added which will bind to the antigen-antibody complex. Any unbound reagent is removed with a further wash step and the specimen is then visualised with the aid of a fluorescence microscope.

A positive reaction is denoted by bright green fluorescence. Uninfected cells appear dull red due to presence of Evan's blue counterstain in the FITC conjugate.

## Kit Components

### Materials Provided

#### 1. Antigen Control Slides:

SLIDE
-------

5 x 8 well antigen control slides (Cat. No. V4RVPS). Each slide consists of a positive control well infected with each of the following: adenovirus, influenza A & B, parainfluenza 1, 2 & 3 and RSV. The negative control well contains uninfected cells.

#### 2. Monoclonal Antibodies\*\*:

CONTROL	+	A
---------	---	---

1 x 1 ml vial containing monoclonal antibodies against adenovirus.

CONTROL	+	IA
---------	---	----

1 x 1 ml vial containing monoclonal antibodies against influenza A.

CONTROL	+	IB
---------	---	----

1 x 1 ml vial containing monoclonal antibodies against influenza B.

CONTROL	+	P1
---------	---	----

1 x 1 ml vial containing monoclonal antibodies against parainfluenza 1.

CONTROL	+	P2
---------	---	----

1 x 1 ml vial containing monoclonal antibodies against parainfluenza 2.

CONTROL	+	P3
---------	---	----

1 x 1 ml vial containing monoclonal antibodies against parainfluenza 3.

CONTROL	+	R
---------	---	---

1 x 1 ml vial containing monoclonal antibodies against RSV.

#### 3. Antiviral Screening Reagent

SCREEN	+
--------	---

1 x 5 ml of pooled monoclonal antibodies directed against: Adenovirus, influenza A & B, parainfluenza 1, 2 & 3 and RSV.

#### 4. Negative Control:

CONTROL	-	IgG
---------	---	-----

1 x 5 ml of non-immune mouse antibody to be used as a negative control.

#### 5. Anti-Mouse IgG FITC Conjugate\*\*:

CONJ	IgG
------	-----

2 x 5 ml of goat anti-mouse IgG antibody conjugate to fluorescein isothiocyanate (FITC).

#### 6. Mounting Media\*\*:

MM
----

1 x 5 ml of Tris-buffer glycerol containing fluorescence enhancer and sodium azide as a preservative.

#### 7. Wash Buffer Concentrate (PBS)\*\*:

BUF	WASH	CONC
-----	------	------

1 sachet of phosphate buffered saline salts which yields 1 litre when dissolved in distilled water. Store in a clean closed container at room temperature.

8. **Tween 20:**

TWEEN

1 x 5 ml of polyoxyethylene sorbitan monolaurate (Tween 20) and sodium azide (NaN<sub>3</sub>) concentrate to be diluted 1:100 in PBS.

9. **Product Insert:**

INS

Instructions for use.

**\*\* Sodium azide (present in the conjugate, monoclonal antibodies, wash buffer, and mounting fluid) is a potentially biohazardous material. It can react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these materials, flush with large volumes of water to prevent azide build-up.**

***Additional Materials Required***

- Cell culture for isolation of respiratory viruses: each laboratory must maintain viable stocks of cells at appropriate passage-state that will efficiently allow replication of respiratory viruses from processed patient specimens. These cells must be checked periodically for ability to support growth of respiratory viruses. The cell lines most commonly used are MDCK, LLC-MK2, A549, Hep-2 and diploid fibroblasts (W138, HNF, MRC-5).
- Viral transport medium (VTM), which is non-inhibitory to the respiratory viruses and the tissue culture cells, adenovirus, influenza A & B, parainfluenza 1, 2 & 3 or RSV, used for viral stabilisation. Hank's balanced salt solution with antibiotics and protein stabiliser is a suitable medium.

**Note:** animal sera other than precolostral foetal bovine serum should not be used as a protein stabiliser to avoid interference from inherent antibody.

- Tissue culture maintenance medium. RPMI or Eagles Minimal Essential Medium with the appropriate amount of precolostral foetal bovine serum are suitable for maintenance after virus infection.
- Sterile tissue culture tubes, dram vials or multi-well plates.
- Sterile swabs (cotton or dracon) which are non-inhibitory to respiratory viruses and tissue culture cells.
- Vials for specimen transport and collection.
- Sterile Pasteur pipettes or accurate pipettes and disposable tips.
- Sterile graduated pipettes to deliver 1ml, 5ml and 10ml.
- Fine tipped forceps.
- Sterile glass beads (1-3mm diameter).
- High quality distilled or deionised water.
- Sodium hypochlorite (0.05%) solution.
- Acetone 99.5%

[NOTE: acetone is hygroscopic and should be kept in a tightly sealed container. Acetone contaminated with water may result in a hazy appearance of the substrate in fluorescent assays].

- Glass slides (preferably with 2 and 8 wells for screening and identification) - acetone rinsed and cleaned.
- No. 1 thickness glass coverslips (22 x 50 mm).
- Moist chamber for incubating slides (35-37°C).
- Fluorescence microscope with correct filter combination for FITC (excitation peak = 490nm, emission peak = 520nm).
- Vortex mixer or sonicator.
- Centrifuge.
- Timer
- Wash Bottle
- Incubator with rheostat for temperature regulation (35-37°C)

## Precautions

### Safety

- For *in vitro* diagnostic use only.
- The kit is intended for use by qualified laboratory staff only.
- Conjugate, monoclonal antibodies, wash buffer and mounting medium contain sodium azide as a preservative. Sodium azide may form potentially explosive metal azides with lead / copper plumbing. For disposal, reagent should be flushed with large volumes of water to prevent azide build up.
- The anti-mouse IgG conjugate contains Evans Blue counterstain which is a potential carcinogen. If skin contact occurs, flush with large volumes of water immediately.
- Acetone is extremely flammable and harmful if swallowed or inhaled. Keep away from heat, sparks, or flames. Avoid breathing vapour. Use adequate ventilation.

- Although control slides have been inactivated, they should be handled and disposed of as with other potentially infectious material.
- Dispose of all clinical specimens, infected or potentially infected material in accordance with Good Laboratory Practice. All such materials should be handled and disposed of as though potentially infectious.
- Do not pipette materials by the mouth and never eat or drink at the laboratory workbench.
- Only experienced personnel should attempt tissue culture procedures.
- Residues of chemicals, preparations and kit components are generally considered as hazardous waste. All such materials should be disposed of in accordance with established safety procedures.
- Wear protective clothing, disposable latex gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.

### Procedural

- Do not use kit or individual reagents past their expiry date.
- Alteration of protocol provided may cause erroneous results.
- Do not mix or substitute reagents from different kit lot numbers.
- Do not substitute reagents from other manufacturers.
- Pooling or diluting of conjugates or monoclonal antibodies may cause erroneous results.
- Avoid leaving reagents above 4°C for prolonged periods.
- Do not expose reagents to bright light during storage or incubation.
- Negative Control (normal mouse antibody) should be tested with each cell culture isolate. Fluorescence indicates a non specific reaction and the test is considered invalid.
- Do not allow slides to dry at any time during the staining procedure.

- If staining multiple specimens on the same slide, care should be taken to avoid cross contamination between samples.
- Always use clean, preferably disposable, glassware for all reagent preparation.

### Storage and Stability

- The kit is stable until the expiry date indicated on the outer box label, provided it is stored at between 2-8°C. Do not freeze or expose to elevated temperatures. Discard any remaining reagents after the kit expiration date.
- The conjugate and monoclonal antibodies are stable for 12 months when stored away from light at 4°C.
- PBS should be stored in a clean, closed container at room temperature.
- During incubation, slides should be protected from light and kept in a humid chamber.
- Acetone is hygroscopic and should be kept in a tightly sealed container.

## **Specimen Collection and Preparation**

Correct sample collection, transport, preparation and storage is essential for successful laboratory diagnosis. The possibility of viral isolation is increased when specimens are collected as soon as possible (3-7 days) after the onset of symptoms. Depending on the type of respiratory syndrome experienced by the patient, several different specimens may be collected.

### ***Specimens of Choice***

#### ***Upper Respiratory Tract Infections***

*Colds* - Nasal washes and aspirates (normally give higher titres) but nasopharyngeal swabs are also acceptable.

*Pharyngitis* - Nasopharyngeal swabs or aspirates should be collected (especially if nasal symptoms are prominent). Throat washings or swabs may be preferable if there is prominent pharyngitis.

*Laryngitis* - Nasopharyngeal swabs or aspirates should be collected.

*Croup/Laryngotracheobronchitis* - Nasopharyngeal aspirates are the preferred specimen.

#### ***Lower Respiratory Tract Infections***

*Bronchitis/tracheobronchitis* - Nasopharyngeal aspirates are the preferred specimen.

*Bronchiolitis* - Nasopharyngeal aspirates are the preferred specimen.

*Pertussis like syndrome* - Nasopharyngeal swabs or aspirates should be collected.

*Influenza/syndrome* - Nasopharyngeal swabs or aspirates should be collected.

*Pneumonia* - Nasopharyngeal swabs or throat swabs are acceptable.

## ***Specimen Collection***

*Nasopharyngeal swabs:* A dry swab should be inserted into one or both nostrils into nasopharyngeal area (a separate swab for each nostril may increase specimen volume). Leave the swab to absorb secretions for a few seconds, rotate gently and remove. Place the swab into 1-2 ml viral transport medium, break off the shaft and seal container tightly.

*Throat swabs:* Moisten a swab in sterile PBS or VTM and rub rigorously on the tonsils and posterior pharynx. Place the swab into 1-2ml of VTM, break off the shaft and seal container tightly.

*Throat wash:* The patient should gargle with 3-5 ml of sterile PBS, which is then collected into a sterile vial containing 1-2ml of VTM and seal tightly.

*Nasal wash:* 2-3ml sterile PBS is flushed into each nostril while the patient holds their head tilted back. The head is then tipped forward and the fluid

collected into a sterile vial containing VTM and seal tightly.

*Nasopharyngeal wash:* 3-7ml of sterile PBS is aspirated into a soft rubber bulb. Place the patient on their side and gently close one nostril with finger pressure. Use the point of the bulb to completely occlude the other side and then squeeze PBS into the nose and aspirate quickly. Secretions are then expelled into a sterile vial containing 1-2ml of VTM and tightly sealed.

*Nasopharyngeal aspirates:* A no. 8 French soft plastic feeding tube is connected through a valve-containing-trap to an electric suction pump apparatus. Introduce the sterile catheter through the nares to the back of the nose. Apply suction intermittently while slowly withdrawing the catheter. The process may be repeated once in each nostril so 0.2-0.8ml of secretion is obtained. The secretion should then be transferred to 1-2ml of VTM and the container sealed tightly.

### ***Specimen Transport and Storage***

All specimens should be transported to the laboratory on ice immediately after collection. Most respiratory viruses are unstable and sensitive to repeated freezing and thawing. Therefore, specimens should be held at 2-8°C prior to inoculation. If specimens cannot be processed within 72 hours, freezing at -70°C or below is recommended. A quick freeze in an acetone/dry ice bath will help maintain viral infectivity. However, freezing specimens greatly reduces the chance of isolation.

### ***Specimen Preparation***

- In order to dislodge cells trapped within the swab fibres, swabs should be rotated while immersed in the transport medium. Swabs should then be discarded into sodium hypochlorite solution.
- For increased cell recovery, add a few sterile glass beads to the specimen and vortex or sonicate (8-12kc/sec) for up to 1 minute.

- Centrifuge the specimen at 2000xg for 10 minutes to remove bacterial contaminants and debris. The supernatant should be used as the inoculum material.

### ***Processing Specimens for Direct Examination***

If specimens are to be used for both direct detection and cell culture confirmation half of the cells should be removed by centrifugation at 300 to 500xg and used with the following procedure. Supernatant and remaining cells from the specimen should be used under the procedure for "Cell Culture Isolation and Fixation," (see page 14).

1. Remove specimen material from original container and place in a 10-15ml sterile centrifuge tube.
2. Add 4ml of sterile PBS to the specimen and mix.
3. Centrifuge specimen at 300-500xg for 10 minutes at 2-8°C.

4. The supernatant is collected to be used for viral isolation procedures.

[NOTE: Refer to "Cell Culture Isolation and Fixation" section for complete instructions. (see page 14)].

5. The cell pellet is washed by adding 4-8ml of PBS and gently resuspending cells.
6. Centrifuge at 300-500xg for 10 minutes at 2-8°C.
7. If mucous is present, it will form a hazy layer above the cell pellet. Carefully remove supernatant and mucous with a sterile Pasteur pipette. If mucous is still present, repeat steps 6-8 until mucous is removed.
8. Resuspend cells in 0.1-0.2ml PBS to make a cloudy suspension. Excessively thick suspensions are difficult to read and do not make high quality slides. Cell suspensions that do not contain enough cells result in a loss of sensitivity. Smears should contain at least two epithelial cells per 250x field.

9. Place one drop of cell suspension onto desired number of wells on pre-cleaned slides. Allow slide to air dry.
10. Fix slides in chilled (2-8°C) acetone for 10 minutes. Do not allow acetone to become contaminated with water and salts. This can produce hazy staining.
11. Allow slides to air dry after fixation. The slides should be stained as soon as possible. If storage is necessary, place slides in desiccated container at -20°C. Slides kept at -20°C will keep for up to one year. Slides must be stored in air-tight containers to prevent moisture penetration.
12. Proceed to "Specimen Screening" section, (see page 15).

### ***Cell Culture Isolation and Fixation***

1. Examine the cell cultures immediately prior to inoculation to ensure good cell viability and morphology.

2. Remove old medium from cell culture using a sterile pipette and refeed the monolayer with culture maintenance medium. Add at least 2ml of fresh medium to glass culture tubes (16x125mm) or add at least 1ml of medium to dram vials.

3. Add 0.2-0.5ml of the specimen to each culture vessel and incubate for 60 minutes at 35-37°C. It is recommended that specimens are inoculated in duplicate.

[NOTE: viral adsorption of the specimen and/or low speed centrifugation may enhance isolation of some viruses and CPE may appear as early as 48 hours].

4. To ensure infection sensitivity and appropriate CPE, a representative (commercial) virus strain should be used to inoculate the cell line. After the inoculation period, cover the cells with fresh maintenance media and incubate in a stationary position at 35-37°C.

5. The culture medium should be renewed every 3 to 4 days to enhance CPE formation.

6. Examine the cells on a daily basis for CPE.

7. When CPE is observed, fix the cells (as follows) in preparation for staining confirmation.

8. Aspirate the culture medium using a sterile pipette and store in a sterile tube. (Subculturing can be carried out from this medium if the cell culture is destroyed during staining).

9. Rinse the cells gently 3 times with 1-2ml volumes of Hank's balanced salt solution. Discard rinses into sodium hypochlorite solution.

10. Add one tenth the culture volume of trypsin (0.05%) Sodium – EDTA (0.53mM) and leave to stand for 30 seconds.

11. Gently tap the culture vessel to loosen the cells.

12. Add sufficient fresh medium to return to the original cell culture volume.

13. Centrifuge the cell suspension at 200-500xg for 7-10minutes.

14. Resuspend the cell pellet in several drops of sterile PBS to give a relatively heavy suspension (approx.  $2 \times 10^6$  cells/ml).

15. Spot the cell suspension onto at least one screening (2 wells) and one identification slide (8 wells) and dry in a convection oven (30-35°C) or air dry rapidly at room temperature.

[NOTE: 5 mm diameter wells are recommended]

16. Fix the slides in chilled acetone (2-4°C) for 10 minutes and air dry completely.

17. Store unused slides with desiccant at -20°C.

### Assay Procedure

#### Reagent Preparation

Preparation of Wash Buffer (PBST):

1. Dissolve the contents of the PBS sachet in 990ml of distilled water.
2. To 495ml of this solution, add 5ml of Tween 20/sodium azide (100x concentrate) and mix thoroughly.

[Note: the remaining PBS solution should be discarded].

3. Transfer to a clean labelled storage container and seal cap tightly.

All remaining reagents are **ready-to use** and are at working dilution.

BIOTRIN Respirator y Virus Panel	NEG	AD	IA	IB
	P1	P2	P3	RSV

Figure 1: Sample antigen control slide

#### Specimen Screening

1. Remove specimen slides (2 wells) and one antigen control slide from storage and allow them to equilibrate to room temperature.

2. Add 1 drop (40µl) of anti-viral screening reagent to one well of the specimen screening slide and each antigen control well of the antigen control slide. [NOTE: the anti-viral screening reagent is used to determine the presence of adenovirus, influenza A, influenza B, parainfluenza 1, parainfluenza 2, parainfluenza 3 or RSV].
3. Add 1 drop (40µl) of normal mouse antibody to a second well of the specimen screening slide and to the negative well of the antigen control slide. [NOTE: this negative control allows detection of any non specific fluorescence].
4. Incubate the specimen and control slides for 30 minutes at 37°C in a humid chamber.
5. Rinse slides thoroughly with Wash Buffer for 10-15 seconds. Shake off any excess wash buffer from the slides and dry between the wells.

6. Add 1 drop (40µl) of anti-mouse IgG FITC conjugate to each well on the specimen and control slides.
7. Incubate the specimen and control slides for 30 minutes at 37°C in a humid chamber.
8. Rinse slides thoroughly with Wash Buffer solution for 10-15 seconds. Shake off any excess wash buffer from the slides and dry between the wells.
9. Mount the slides with no. 1 coverslips and mounting media. Remove excess mounting media from the slides.
10. Examine the slides under a fluorescence microscope at 100-200x. Detailed examination can then be carried out at 400x. [NOTE: the control slide should be examined first for appropriate staining].
11. If screening slides are negative, they should be reported as "no virus observed". Proceed to step 12 for all specimens that are positive.

### ***Virus Identification***

1. Remove the identification slides (8 wells) corresponding to each positive specimen and an antigen control slide from storage. Allow slides to equilibrate to room temperature.
2. Add 1 drop (40µl) of each monoclonal antibody to each well on the specimen identification slide.
3. Add 1 drop (40µl) of each monoclonal antibody to the corresponding well on the antigen control slide.
4. Add 1 drop (40µl) of Negative Control to a separate well on the specimen identification slide and to the negative well of the antigen control slide.
5. Incubate the specimen and control slides for 30 minutes at 37°C in a humid chamber.
6. Rinse slides thoroughly with wash buffer for 10-15 seconds. Shake off any excess wash buffer from the slides and dry between wells.

7. Add 1 drop (40µl) of anti-mouse IgG FITC conjugate to each well on the specimen and control slides.
8. Incubate the specimen and control slides for 30 minutes at 37°C in a humid chamber.
9. Rinse slides thoroughly with wash buffer for 10-15 seconds. Shake off any excess wash buffer from the slides and dry between the wells.
10. Mount the slides with no. 1 coverslips and glycerol mounting media. Remove excess mounting media from the slides.  
[NOTE: slides may be stored in the dark, in an air tight box at 2-8°C for 24 hours without significant loss of fluorescence. However, immediate observation after staining is recommended. For longer storage, slides should be held in an air tight box at -20°C or lower].

## Interpretation of Results

Typical staining patterns exhibited by the infecting virus are described below. The type of staining seen on the slide will be dependent on the virus/viruses present and reflects its growth pattern.

### *Adenovirus*

Bright green fluorescence is nuclear, cytoplasmic, or both. Nuclear staining is uniformly bright with little definition. Cytoplasmic staining is often punctuate.

### *Influenza A and B*

Bright green fluorescence is nuclear, cytoplasmic, or both. Nuclear staining is uniformly bright with little definition. Cytoplasmic staining is often punctuate with large inclusions.

### *Parainfluenza 1, 2, and 3*

Bright green fluorescence is confined to the cytoplasm. Cytoplasmic staining is punctuate with irregular inclusions.

### *Respiratory Syncytial Virus*

Bright green fluorescence is seen in the cytoplasm and associated with syncytia. Cytoplasmic staining is punctuate with small inclusions.

### *Negative Cells*

Cells stain with a dull red in the cytoplasm and a deep maroon almost black in the nucleus. Staining of these cells is due to the presence of Evans Blue counterstain.

It is recommended that negative cells are examined first to establish if there is any background fluorescence present.

A sample is considered positive for a virus if at least 2 or more bright green fluorescent cells are observed per field at 400x magnification.

The positive control should show cells exhibiting fluorescence of the nucleus, cytoplasm or both, depending on the stage of growth. The negative control should show a dull red staining due to the presence of Evans Blue counterstain.

[NOTE: Fluorescent staining of cell fragments should be ignored; this may occur due to the trapping of the conjugate in such debris. If the screening slide shows a positive result and the identification slide shows no positive result, the test should be repeated to confirm results].

## Expected Values

Isolation and identification results will vary depending on geographic location, age of population, time of year (some viruses are seasonal), socio-economic status, specimen handling and the type of detection system used.

For the period September to July (1991-1992) the following detection rates were obtained: \*

Respiratory Virus	Detection rates (sample no.)	Detection rates (%)
Adenovirus	43/646	6.7%
Influenza A	11/646	1.7%
Influenza B	4/646	0.6%
Parainfluenza 1	11/646	1.7%
Parainfluenza 2	5/646	0.8%
Parainfluenza 3	28/646	4.3%
RSV	66/646	10.2%

\*Settings for the evaluations included hospitals, children's hospitals, universities, and reference laboratories in the USA.

## Limitations of Use

- Positive and negative controls should be tested with each specimen and must give the appropriate staining for the test to be valid. [NOTE: additional control slides are available - Cat No: V4RVPS (5 per pack)].
- A negative result does not preclude the presence of adenovirus, influenza A & B, parainfluenza 1, 2 & 3, or RSV. Results should be interpreted in conjunction with patient related clinical information and other diagnostic procedures. Failure to detect virus may be the result of factors such as improper specimen collection or handling, improper culture techniques or the use of an inappropriate cell line or temperature during isolation. All negative results should be reported as "no virus observed".

- The monoclonal antibodies in this kit are group specific for adenovirus and RSV; therefore, they cannot be used for type differentiation.

*Staphylococcus aureus* contaminated samples may exhibit dull yellow-green fluorescence due to the presence of large amounts of protein A. The fluorescence is a result of non-specific binding of protein A to antibody F<sub>c</sub> fragments.

## Performance Characteristics

- Clinical Comparison of Cell Culture Confirmation*

The following results were obtained when cell culture confirmation was carried out on respiratory samples (a total of 646 specimens were tested):

MONOCLONAL ANTIBODIES:	Adeno	RSV	Inf. A	Inf. B	Para 1	Para 2	Para 3
Total Positive	43	66	11	4	11	5	28
Sensitivity %	97.7	100	100	100	100	100	100
Specificity %	100	100	100	100	100	100	100
Predictive Value (+) %	100	100	100	100	100	100	100
Predictive Value (-) %	99.8	100	100	100	100	100	100

The following results were obtained when cell culture confirmation was carried out on frozen respiratory samples (a total of 151 specimens were tested):

MONOCLONAL ANTIBODIES:	Adeno	RSV	Inf. A	Inf. B	Para 1	Para 2	Para 3
Total Positive	7	3	30	36	33	5	16
Sensitivity %	100	100	100	100	100	100	100
Specificity %	100	100	100	100	100	100	100
Predictive Value (+) %	100	100	100	100	100	100	100
Predictive Value (-) %	100	100	100	100	100	100	100

- Cross Reactivity

The following results were obtained from viral and bacterial detection studies with each monoclonal antibody.

MONOCLONAL ANTIBODIES:	Adeno	RSV	Inf. A	Inf. B	Para 1	Para 2	Para 3
Adenoviruses	+	-	-	-	-	-	-
Coronaviruses	-	-	-	-	-	-	-
Coxsackievirus	-	-	-	-	-	-	-
Echoviruses	-	-	-	-	-	-	-
Herpes viruses	-	-	-	-	-	-	-
Influenza virus Type A	-	-	+	-	-	-	-
Influenza virus Type B	-	-	-	+	-	-	-
Measles virus	-	-	-	-	-	-	-
Mumps virus	-	-	-	-	-	-	-
Parainfluenza Type 1	-	-	-	-	+	-	-
Parainfluenza Type 2	-	-	-	-	-	+	-
Parainfluenza Type 3	-	-	-	-	-	-	+
Rhinoviruses	-	-	-	-	-	-	-
Simian virus	-	-	-	-	-	-	-

RSV	-	+	-	-	-	-	-
Acholeplasma laidlawii	-	-	-	-	-	-	-
Bordetella bronchiseptica	-	-	-	-	-	-	-
Bordetella pertussis	-	-	-	-	-	-	-
Bordetella parapertussis	-	-	-	-	-	-	-
Branhamella catarrhalis	-	-	-	-	-	-	-
Chlamydia trachomatis	-	-	-	-	-	-	-
Corynebacterium diphtheriae	-	-	-	-	-	-	-
Legionella micdadei	-	-	-	-	-	-	-
Legionella pneumophila	-	-	-	-	-	-	-
Mycobacterium avium	-	-	-	-	-	-	-
Mycobacterium intracellulare	-	-	-	-	-	-	-
Mycobacterium tuberculosis	-	-	-	-	-	-	-
Mycoplasma fermentans	-	-	-	-	-	-	-
Mycoplasma hominis	-	-	-	-	-	-	-
Mycoplasma orale	-	-	-	-	-	-	-
Mycoplasma pneumoniae	-	-	-	-	-	-	-
Neisseria meningitidis	-	-	-	-	-	-	-
Ureaplasma urealyticum	-	-	-	-	-	-	-

Host Cell Controls	All monoclonal antibodies tested negative against the following cell controls
Adenovirus	pHEK, HEp-2, NCI-H292, RD, HLF, HeLa, & A549 cell controls
Polyomaviruses	PMK cell control
Enteroviruses	RD & HLF cell controls
Paramyxoviruses	HEp-2, Vero, & PMK cell controls

## Summary of Respiratory Virus Panel Procedure


Important Note: Please read the entire product instruction leaflet before starting the assay. This summary is for quick reference only.

### Specimen Screening

Add 1 drop (40µl) of anti-viral screening reagent to one well of the specimen screening slide and each antigen control well of the antigen control slide

Add 1 drop (40µl) of normal mouse antibody to a second well of the specimen screening slide and to the negative well of the antigen control slide

Incubate for 30 minutes @ 37°C in a humid chamber

## Summary of Respiratory Virus Panel Procedure


Important Note: Please read the entire product instruction leaflet before starting the assay. This summary is for quick reference only.

### Virus Identification

Add 1 drop (40µl) of each monoclonal antibody to each well on the specimen identification slide.

Add 1 drop (40µl) of each monoclonal antibody to the corresponding well on the antigen control slide.

Add 1 drop (40µl) of normal mouse antibody to a separate well on the specimen identification slide and to the negative well of the antigen control slide.

## Interpretation of Symbols

*In-vitro* diagnostic medical device 

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1. **Wadell, G.** (1987) Adenoviruses. In principles and practice of clinical virology (Eds. Zuckerman et al) John Wiley and Sons Ltd.
2. **Horowitz, M.S.** (1985) Adenovirus diseases. In Virology 9. (Eds. Fields B.N. et al.) Raven Press New York pp 477-495.
3. **Hierholzer, J.** Adenoviruses. Ch. 86 In : Manual of Clinical Microbiology. 5th Ed. (Eds Balows et al). ASM Washington D.C. pp 896.
4. **Bell J. et al** (1955) Pharyngoconjunctival fever. Epidemiological studies of a recently recognised disease entity. JAMA 175, pp1083-1092.
5. **Kasel, J** (1980) Adenoviruses Ch. 9 In : Viral Rickettsial and Chlamydial Infections. American Publ. Health Assos. (Eds. Lennett, E. and Schmidt, N.)
6. **Cooney, MK.** (1985) Adenoviruses. In Manual of Clinical Microbiology, 4th Ed. (Eds. Lennett, E. et al) American Society for Microbiology, Washington.
7. **Harmon M.W. and Kendal A.P.** Influenza Viruses. Ch. 81 In : Manual of Clinical Microbiology. 5th Ed. (Eds Balows et al). ASM Washington pp868.
8. **Glezen WP, Paredes A, Taber LH.** Influenza in children. Relationship to other respiratory agents. JAMA 1980; 243: 1345-1349.

9. **Frank, A.L et al** (1979) Comparison of different tissue cultures for isolation and quantitation of influenza and parainfluenza viruses. *J. Clin. Microbiol.* 10
10. **McIntosh, K and Chancock, R.M.** (1990) Parainfluenza viruses. Ch 35 In *Fields Virology Vol. 1* (Eds Fields B.N and Knipes D.M). Raven Press New York, p 963.
11. **McIntosh, K and Chancock, R.M.** (1990) Respiratory syncytial virus. Ch 38 In *Fields Virology Vol. 1* (Eds Fields B.N and Knipes D.M). Raven Press New York, p 1045.
12. **Zinsser Microbiology.** (Eds. Joklik W.K. et al) 19th Edn. Appleton & Lange.
13. **Gardner, P.S. and McQuillin** (1980) *Rapid Virus Diagnosis- applications of immunofluorescence.* 2nd Ed. Butterworth and Co. Ltd. London.
14. **Ahulwalia, G et al.** (1987) comparison of nasopharyngeal aspirate and nasopharyngeal swab specimens for respiratory syncytial virus diagnosis by cell culture, indirect immunofluorescence assay, and enzyme-linked immunosorbent assay. *J. Clin. Microbiol* 25:763-767.
15. **Chernesky, M.A et al** (1982) *Laboratory diagnosis of viral infections.* Coordinating Ed. Drew E.L. American Society for Microbiology, Washington.
16. **Hall, C.B.** (1985) Acute laryngotracheobronchitis and croup p. 360-364 in *Principles and practice of infectious diseases.* (Eds. Mandell G et al.). Wiley and Sons, Inc. New York.